“酷柏” 佰美視雙週拋軟式隱形眼鏡
“CooperVision” Biomedics 55 (Ocufilecon D) Soft Contact Lens

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**BIOMEDICS® 55 (ocufilcon D)
Soft (Hydrophilic) Contact Lenses
Package Insert**

IMPORTANT: Please read carefully and keep this information for future use. This package insert is intended for the eye-care practitioner but should be made available to patients upon request. The eye-care practitioner should provide the patient with the appropriate patient instructions and/or patient information booklet which pertains to the patient’s prescribed lenses.

Symbol	Description
	CAUTION: Federal (U.S.A.) Law restricts this device to sale by, or on the order of a licensed practitioner.
	See Instruction Leaflet
	Use by date (expiration date)
	Batch code
	Sterile using steam heat

INDICATIONS (USES):

BIOMEDICS® 55 UV Asphere (Ocufilcon D) Soft contact lenses is indicated for the correction of visual acuity in persons with non-aphakic, non-diseased eyes which manifest myopia (nearsighted), hyperopia (farsighted), or astigmatism correction lower than –2.00 diopters that does not interfere with visual acuity . It is a daily wear and 2 weeks replacement.

DESCRIPTION:

The lens material for BIOMEDICS 55 UV ASPHERE (ocufilcon D) Soft Contact Lens is a random copolymer of 2-hydroxyethylmethacrylate and methacrylic acid. The BIOMEDICS 55 UV ASPHERE (ocufilcon D) Soft Contact Lenses contain visibility blue tint from edge to edge using (VAT Blue 6) which is added in monomer prior to polymerization to make the lens more visible for handling. The BIOMEDICS 55 UV ASPHERE (ocufilcon D) Soft Contact Lenses contain a benzophenone UV absorbing monomer which is used to block UV radiation. The BIOMEDICS 55 UV ASPHERE (ocufilcon D) Soft Contact Lenses help protect against transmission of harmful UV radiation to the cornea and into the eye.

BIOMEDICS 55 UV ASPHERE (ocufilcon D) Soft contact lenses is hemispherical shells of the following dimensions:

Diameter:	12.0 mm to 18.0mm
Base Curve:	6.50 mm to 12.80 mm
Center Thickness:	0.025 mm to 0.40 mm,
Powers:	Plano to – 10.00D


The physical parameters of the lens are as follows:

Refractive Index:	1.41
Surface Characteristics:	Hydrophilic
Water Content:	55.0%
Oxygen permeability (Dk)	19.6 X 10 ⁻¹¹ (cm ² /sec) (ml O ₂ /ml x mm Hg at 35°C
Light Transmittance:	97.0%

Design

The BIOMEDICS 55 UV ASPHERE (ocufilcon D) Soft Contact Lens is hemispherical flexible shells which cover the cornea and may cover apportion of the adjacent sclera.

This is the transmittance profile of –8.00D ocufilcon D hydrophilic contact lens with UV blocker versus a human cornea and a human crystalline lens

1. Lerman, S., *Radiant Energy and the Eye* (New York: MacMillan, 1980), p. 58, Figure 2-21. Transmittance profile of the human cornea of a 24-year-old person.
2. Waxler, M., and V. M. Hitchens, *Optical Radiation and Visual Health* (Boca Raton: CRC Press, 1986), p. 19, Figure 5. Transmittance profile for the human crystalline lens of a 25-year-old person.

WARNING: UV-absorbing contact lenses are NOT substitutes for protective UV-absorbing eyewear, such as UV-absorbing goggles or sunglasses, because they do not completely cover the eye and surrounding area. Patients should continue to use UV-absorbing eyewear as directed.

ACTIONS:

In its hydrated state, the **BIOMEDICS® 55** UV Asphere (Ocufilcon D) Soft contact lens, when placed on the human cornea, acts as a corrective refracting medium to focus light rays on the retina.

The visibility-tinted **BIOMEDICS® 55** UV Asphere (Ocufilcon D) Soft contact lens, allows the lens to become readily visible to the wearer when it is not on the eye.

The thinnest **BIOMEDICS® 55** UV Asphere (Ocufilcon D) Soft contact lens,(-2.00 to -10.00 D) blocks 70% of UVA radiation and 96% UVB radiation average across the spectrum. The radiation blockage of the Biomedics UV lenses will increase for thicker lenses (Please refer to accompanying transmittance curve graph).

NOTE: Long-term exposure to ultraviolet radiation is one of the risk factors associated with cataracts. Exposure is based on a number of factors, such as environmental conditions (altitude, geography, cloud cover) and personal factors (extent and nature of the outdoor activities). UV-absorbing contact lenses help provide protection against harmful UV radiation. However, clinical studies have not been done to demonstrate that wearing UV-absorbing contact lenses reduces the risk of developing cataracts or other eye disorders. Consult their eye-care practitioners for more information.

CONTRAINDICATIONS (REASONS NOT TO USE):

Patients **SHOULD NOT USE BIOMEDICS® 55** UV Asphere (Ocufilcon D) Soft contact lens when any of the following conditions exist:

- Acute and subacute inflammation or infection in the anterior chamber of the eye.
- Any eye disease, injury, or abnormality affecting the cornea, conjunctiva, or eyelids.
- Severe Insufficiency of lachrymal secretion (dry eyes).
- Corneal hypoesthesia (reduced corneal sensitivity).
- Any systemic disease that may affect the eye or be exaggerated by wearing contact lenses.
- Allergic reactions of ocular surfaces or adnexa which may be induced or exaggerated by wearing contact lenses or by using contact lens solutions.
- Allergy to any ingredient, such as mercury or thimerosal or chlohexidine, in a solution which is to be used to care for the lens.

- Use of a medication that is contraindicated, including eye medications.
- Patient history (i.) of recurring eye or eyelid infections, including sties; (ii.) of adverse effects associated with contact lens wear; or (iii.) of intolerance or abnormal ocular response to contact lens wear.
- History of patient non-compliance with contact lens care and disinfection regimens, wearing restrictions, wearing schedule, or follow-up visit schedule.
- Patient inability or unwillingness, because of age, infirmity, or other mental or physical conditions, or because of an adverse working or living environment, to understand or comply with any warnings, precautions, restrictions, or directions. Additionally, patients who require only vision correction and (i.) who would not, or could not, adhere to a recommended care system for lenses; or (ii.) who are unable to place or remove lenses should not be provided with them.

PRECAUTIONS:

- Contact lens wear may not be suitable for those in certain occupations, or, in other instances, such persons may require protection equipment.

- In order to minimize the likelihood of lens contamination or of physical trauma to the cornea, lens-wearing patients should avoid environmental fumes, smoke, dust, vapors, and windy conditions.

- Recommend to use of sterile lens-care solutions. If a particular patient is allergic to preservatives, that patient should use sterile non-preserved solutions and should discard such solutions after the time specified in their label directions.

- Eye injury from irritation or infection and damage to contact lenses may result from lens contamination. Patients should take care to prevent cosmetics, lotions, soaps, creams, hair sprays, or deodorants from coming into contact with their lenses.

- Tweezers or other tools should not be used by patients to remove lenses from lens containers; rather, the contents of a lens container should be poured into the hand.

- Practitioners should instruct their patients as to the proper manner to promptly remove their lenses, and patients should be able to demonstrate the ability to do so.

- Fluorescein should not be used while **BIOMEDICS® 55** UV Asphere (Ocufilcon D) Soft contact lens is on the patient’s eye. The lenses absorb this dye and become discolored. In the event fluorescein does come in contact with the lenses while they are on the patient’s eye, the eyes should be flushed thoroughly with a sterile saline solution recommended for in-eye use, and new lenses should be inserted only after at least one hour.

- After removal of the lenses from the lens case, to prevent contamination and to help avoid serious eye injury, the patient should always empty and rinse the lens case with fresh rinsing solution and allow it to air-dry between each lens disinfection cycle.

- Oral contraceptive users could develop visual changes or changes in lens tolerance when using contact lenses. Patients should be cautioned accordingly.

- Contact lens-wearing patients should be instructed to inform their physicians that they wear contact lenses; further, patients’ physicians should consult their eye-care practitioners before using any medication in the eye.

- Certain medications (such as antihistamines, decongestants, diuretics, muscle relaxants, tranquilizers, and medications for motion sickness) may cause dryness of the eye, increased lens awareness, or blurred vision. Should these conditions exist, proper remedial measures should be prescribed. Depending on the severity, such measures could include the use of lubricating drops that are indicated for use with soft contact lenses or the temporary discontinuance of contact lens wear while such medications are being used.

- As with any contact lens, follow-up visits are necessary to better ensure the continuing health of the patient’s eyes. The patient should be instructed as to a recommended follow-up schedule.

- Practitioners should caution their patients wearing **BIOMEDICS® 55** UV Asphere (Ocufilcon D) Soft contact lens on a daily-wear schedule to remove their lenses before sleeping.

ADVERSE EFFECTS:

Patients should be informed that the following problems may occur when they wear contact lenses:

- The eye may be painful or may burn, sting, or itch.
- There may be less comfort than when the lens was first placed on eye.
- There may be a feeling of something in the eye (foreign body, scratched area).
- There may be the potential for some temporary impairment due to peripheral infiltrates, peripheral corneal ulcers, and corneal erosion. There may be the potential for other physiological observations, such as local or generalized edema, corneal neovascularization, corneal staining, infection, tarsal abnormalities, iritis, and conjunctivitis, some of which are clinically acceptable in low amounts.
- There may be excessive watering (tearing), unusual eye secretions, or redness of the eye.
- Poor visual acuity; rainbows or halos around objects; photophobia; or a feeling of dryness (dry eyes) may also occur if the lenses are worn continuously or for too long a time.
-

If a patient reports any problems, including, but not limited to, the foregoing, he or she should be instructed to remove his or her lenses immediately. Then:

- If the discomfort or problem stops, the patient should look closely at the lens.
- If the lens has dirt, an eyelash, or other foreign body on it, the patient should be instructed as follows:
 - If the patient is in the Disposable Wear Program, he or she should dispose of the lens and replace it with a fresh, new lens.
 - If the patient is in the Scheduled Replacement Program, and if the lens appears undamaged, he or she may clean, disinfect, and reinsert the lens.
- If the lens is or appears in any way damaged, the patient SHOULD NOT put the lens back on the eye. The patient should discard the lens and insert a fresh, new lens on the eye.
- If the patient’s problem continues, the patient SHOULD NOT put the lens back on the eye; but rather he or she should immediately contact his or her eye-care practitioner or a physician, who must determine the need for examination, treatment, or referral without delay.
- The patient should be advised that when any of the aforementioned symptoms occur, a serious condition such as infection, corneal ulcer, corneal neovascularization, or iritis may be present and may progress rapidly. The patient should seek immediate professional identification of the problem and prompt treatment to avoid serious eye damage. Less serious reactions, such as abrasions, epithelial staining, and bacterial conjunctivitis, should be treated appropriately to avoid complications.

WEARING SCHEDULES:

It is recommended that a contact lens-wearing patient see his or her eye-care practitioner twice each year or, if so directed, more frequently. The practitioner should determine the appropriate wearing schedule and replacement schedule, which he or she should provide to the patient.

Daily Wear: There may be a tendency for the daily-wear patient to overwear the lenses initially. Therefore, practitioners should stress to these patients the importance of adhering to a proper initial daily wearing schedule. The practitioner should determine the appropriate wearing schedule and replacement schedule, which he or she should provide to the patient.

LENS CARE DIRECTIONS:

Eye-care practitioners should provide their patients with appropriate and adequate instructions and warnings for lens care and handling, and practitioners should recommend appropriate and adequate procedures and products for each individual patient in accordance with the particular lens-wearing schedule and care system selected by the practitioner, the specific instructions for such products, and the particular characteristics of the patient.

For patients in the Disposable Wear Program: Eye-care practitioners should review with patients that no cleaning or disinfection is needed with disposable lenses. Patients should always dispose of lenses when they are removed and have replacement lenses or spectacles available.

For complete information concerning emergency lens care, refer to the Patient Information Booklet for patients in the Disposable Wear Program. Emergency lens care does not apply to lenses worn on a daily-wear basis.

For patients in the Scheduled Replacement Program: For complete information concerning the care, cleaning, and disinfecting of **BIOMEDICS® 55** UV Asphere (Ocufilcon D) Soft contact lens, patients should refer to the Patient Information Booklet for the Scheduled Replacement Program.

CARE FOR A DEHYDRATED LENS:

If a soft contact lens is off the eye and is exposed to air for a significant period of time, it may become dry and brittle and need to be re-hydrated. If the lens is adhering to a surface such as a counter top, apply sterile saline before handling the lens.

Eye-care practitioners should review the following information on re-hydrating the lens with the patient:

- Handle the lens carefully.
- Place the lens in a storage case and soak the lens in a recommended rinsing and storing solution for at least an hour until it returns to a soft state.
- Clean and disinfect and re-hydrate lens using a recommended lens-care system.
- If the lens does not become soft after soaking, the lens should not be used until it is examined by the eye-care practitioner.

CARE FOR A STICKING LENS:

If the lens sticks (or stops moving), the patient should be instructed to apply several drops of the recommended lubricating or rewetting solution directly to the eye and wait until the lens begins to move freely on the eye before attempting to remove it. If non-movement of the lens continues after several minutes, the patient should **immediately** consult with his or her eye-care practitioner immediately.

PRACTITIONER FITTING SETS:

All lenses which have been opened must be discarded after each fitting.

HOW SUPPLIED:

Each **BIOMEDICS® 55** UV Asphere (Ocufilcon D) Soft contact lens is supplied sterile in a container with a normal buffered saline solution. Several containers are packaged in a multi-pack arrangement, each of which is marked with the manufacturing lot number of the lens, the dioptric power, the base curve or series, the diameter, the axis, the cylinder and the expiration date.

Package : 1 lens/box, 2 lenses/box, 3 lenses/box, 6 lenses/box

REPORTING OF ADVERSE REACTIONS:

All serious adverse experiences and adverse reactions observed in patients wearing **BIOMEDICS® 55** UV Asphere (Ocufilcon D) Soft contact lens or experienced with the lenses should be reported

CooperVision Caribbean Corporation
500 Road 584, Lot 7,
Amuelas Industrial Park
Juana Diaz, PR 00795, U.S.A

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